

PURINOX HME[®]

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Figures

Disclaimer

Purinox HME B.V. makes no express or implied warranty to the purchaser hereunder with regard to the lifetime of the delivered product, which may vary according to individual use and biological conditions. In addition, Purinox HME B.V. makes no warranty as to the merchantability or fitness of the product for any particular purpose.

Content

EN - English..... 4

NL - Nederlands..... 8

EN - English

1 Descriptive information

1.1 Intended use

The PURINOX HME® Sample is a 22 mm passive heat and moisture exchanger for patients breathing through a tracheostoma. The PURINOX HME® Sample is intended for single use and for adults only.

1.2 Contra-indications

The device should not be used by patients with impaired cognitive or physical abilities. The device should not be used by patients who are unable to attach or remove the device themselves, or who do not have sufficient knowledge of how to use the device, or who are not cognitively able to understand the risks associated with its use.

The device should not be used by patients with low resting tidal volume, as the additional dead volume may cause retention of CO₂ (carbon dioxide).

1.3 Device description

The PURINOX HME® Sample can be used for a maximum of 24 hours. The PURINOX HME® Sample may be rinsed with water and dried, if necessary, during the 24 hours.

The PURINOX HME® Sample is affixed over the tracheostoma on an tracheostomy adhesive with a 22 mm coupling ring or a tracheostomy tube with a 22 mm connector. A PURINOX HME® cover can be clicked onto the PURINOX HME® Sample to adjust the appearance.

1.4 Technical data

| | | | |
|---|-----------------|-----------------|-----------------|
| Internal volume (ml)* | 4.7 | | |
| Moisture loss at V_T = 1000 ml (mg/l)* | 18.2 | | |
| Pressure drop (hPa)* | 30 l/min | 60 l/min | 90 l/min |
| Initial | 0.4 | 1.0 | 2.0 |
| After 24 hours | 0.4 | 1.3 | 2.3 |

* According to ISO 9360

Material PURINOX HME® Sample: stainless steel 316L

Material PURINOX HME® cover: high density polyethylene

1.5 WARNINGS

- Do not apply unintended force on the PURINOX HME® Sample. Accidental closure of the device may obstruct breathing.
- Do not apply excessive force while attaching the PURINOX HME® Sample or while it is attached. Excessive force may result in the device being pushed through the coupling ring of the adhesive. This can obstruct breathing and block the airways.
- Do not use the PURINOX HME® Sample with any other attachment devices than advised. The device may go through the coupling ring of the adhesive or get stuck in the attachment device. This can obstruct breathing and block the airways.
- Do not attach a PURINOX HME® Sample that is clogged with mucus as this may obstruct breathing.
- Do not reuse the PURINOX HME® Sample. The device is intended for single use only.
- Do not use incompatible or non-genuine devices as this could result in personal injury or damage to the device.
- Do not use the PURINOX HME® Sample during MRI examination. The device contains metal and will interact with the MRI field.
- Consult your physician if you have ANY adverse reaction to the use of the PURINOX HME® Sample.

1.6 PRECAUTIONS

- Always check if the PURINOX HME® Sample is clean and undamaged prior to use.
- Do not attach a PURINOX HME® Sample that appears slippery due to mucus or other substances, as it will be difficult to attach in a proper position.
- Be careful when attaching the PURINOX HME® Sample at night, as dark surroundings may make it more difficult to attach it in the proper position. It may also be hard to detect if the device is contaminated.
- Do not administer medicated nebulizer treatment over the PURINOX HME® Sample since the medication can become deposited in the device
- Do not use humidifiers or heated humidified oxygen tracheostomy while using the PURINOX HME® Sample. The device will become too wet. If oxygen therapy is required, use non-heated humidified oxygen.

- Do not use the PURINOX HME® Sample while showering or during water activities since water may enter the airways.
- Be careful using the device in extremely hot (e.g. sauna, long sun exposure time) or cold environments. The PURINOX HME® Sample functions best between 36 and 94°F (2 and 34°C).
- Be careful using the PURINOX HME® Sample during extreme activities (e.g., contact sports, motorcycle riding) since it can cause unexpected force onto the device or extreme airflow through the device.

2 Instructions for use

2.1 Operation instruction

2.1.1 PURINOX HME® Sample

Preparation

Inspect if the PURINOX HME® Sample is clean and undamaged. Place the attachment device, according to its instructions for use.

To attach

Hold the PURINOX HME® Sample at the side openings between your thumb and index finger. Press lightly until the device connects to the attachment device. Rotate the device such that the airflow holes can easily be sealed with thumb and index finger. See fig. 1.

To speak

Seal the two side openings of the PURINOX HME® Sample with thumb and index finger at the same time to allow speech. See fig. 2.

To remove

Hold the attachment device. Place the index finger in one side opening and the thumb onto the opposite side opening of the PURINOX HME® Sample. Pull lightly sideways until the device disconnects. See fig. 3.

2.1.2 PURINOX HME® cover

To attach

Click the PURINOX HME® cover onto the device preferably when the device is not affixed over the stoma.

To remove

Removing the PURINOX HME® cover is easiest when the device is not affixed over the stoma. Remove the PURINOX HME® cover by placing your finger behind the indentation and pulling lightly sideways until the cover disconnects from the device.

2.2 Lifetime and disposal

The PURINOX HME® Sample can be used for a maximum of 24 hours. The PURINOX HME® Sample can be returned to the manufacturer for safe disposal. Always follow medical practice and national requirements regarding biohazards when disposing of a used medical device.

The PURINOX HME® cover should be disposed of when showing signs of damage or irremovable contamination. The cover can be disposed of with plastic waste.

3 Reporting

Any serious incident that has occurred in connection with the device must be reported to the manufacturer and the national authority of the country where the user and/or patient resides.

Instruction videos are available on www.youtube.com/@purinoxhme

View the most recent version of these instructions for use on
www.purinoxhme.com

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